

VandaCare Competence



HEALTH MANAGEMENT

Keep track of your health and be in-charged of it.



INSTANT ACCESS

Monitor your health anytime, anywhere, everywhere.



ONE CLICK AWAY

Shape your care and love with just one click.



QUICK MEDICAL REFERENCES

Comprehensive medical references for doctors and medical institutions.


Shape Your Life with



Preventive Care through



VandaCare Solutions



Magicsoft Asia Headquater

21 Bukit Batok Crescent
#09-84 Wcega Tower
Singapore 658065
Tel : +(65) 65157087
Fax : +(65) 65157086
ROC 199103693K
kennytew@magicsoft-asia.com

Indonesia Office

Jl Permata Jingga AA-7
Malang 65143
Indonesia
Tel : +(62) 341415518
vita@magicsoft-asia.com

Malaysia Office

96 Jalan SBC 3
68100 Batu Caves, Selangor
Malaysia
Tel : +(60)-19-6589896
Fax : +(60)-3-61881069
kam@magicsoft-asia.com

"An effective way of recording, accessing and managing

all personal health data (blood pressure, glucose, temperature, BMI,

Oxygen level, Cardio info, Heartbeat, Calorie and etc.) at anytime, anywhere

from smart phone and computer wirelessly over the cloud."


...

Realtime Personal Health Care Management System



TD-3250C CLEVER CHEK 2-IN-1 BLOOD GLUCOSE plus BLOOD PRESSURE ARM CUFF METER

- 2-IN-1 system, with Bluetooth technology
- Large internal memory: 352 measurement results with date and time
- Key tone warning: glucose value > 240 mg/dL (13.3 mmol/L)
- Accuracy: 15mg/dL (0.83mmol/L) with glucose 75mg/L (4.17 mmol/L)
20% with glucose 75mg/L (4.17 mmol/L)
3mmHg or 2% of blood pressure reading



W310B FITNESS SCALE WITH BMI

- Weight scale with BMI, with Bluetooth technology
- Large internal memory: 135 measurement results
- Measurement up to 180kg (397lb) in increments of 0.1kg
- Accuracy: 0.5kg
- Weight displayed unit: kg, lb or stones



IR20B EAR THERMOMETER

- Easy and smart ear thermometer, with Bluetooth technology
- Large internal memory: up to 10 measurement results
- Fever indicator and over range warning
- Accuracy: 0.20C (0.40F) for the range of 36.0C to 39.0C (96.80F to 102.20F)
0.30C (0.50F) for the range of 34.0C to 35.90C (93.20F to 96.60F)
0.30C (0.50F) for the range of 39.10C to 42.20C (102.40F to 108.0F)
- Fast result in 1 second



ECG-102 CARDIO CARE ECG/EKG RECORDER

- Handheld ECG/EKG recorder, with Bluetooth technology
- Small and portable
- Standard sensitivity: 10mm/mv 5%
- Input dynamic range: 2mV
- CMRR: > 95db
- Analog to Digital Conversion: 12 bits
- Sampling frequency: 256 samples/sec



POM-102 HEMOXI CARE PULSE OXIMETER

- Handheld Pulse Oximeter, with Bluetooth technology
- Pulse rate range: 40 ~ 250 bpm
- Measuring: Red (660nm), IR (905nm)
- Accuracy: SPO2 - 70~99% 2 digits
Heart Rate - 40~240 bpm 3 digits
- Output rate: 75 bits/sec
- Sampling Rate: 3 Hz